INFORMATIEBLAD

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INFORMATION SHEET

For Optimal Use of the Lithium-ion Battery in Scoozy C

Scoozy features a lithium-ion battery instead of the standard AGM battery. Here are some handy tips to ensure your Scoozy performs at its best and you encounter no unpleasant surprises.

Key Points Summarized:

- Scoozy calculates the range based on the current energy consumption. This can vary due to factors such as speed, weight, and terrain.
- Keep the battery well-charged for optimal performance. With a charge below 50%, the maximum speed slowly decreases.
- We advice to charge the battery when the range is less than 10 km (and with a double battery pack less than 20 km).
- Store Scoozy preferably in a heated environment to keep the battery at a temperature for optimal performance.

How far can I drive with Scoozy?

The theoretical range of Scoozy with one battery is 50 km and with two batteries 100 km. When starting up, Scoozy calculates the range based on the current energy consumption. This range varies depending on various (environmental) factors such as speed and terrain. The practical range is generally about 85% of the theoretical range.

When should I charge the batteries?

When the range with 1 battery pack falls below 10 km and with a double battery pack below 20 km, the range prediction is no longer reliable, and it is advisable to charge the battery.

The advice is to keep Scoozy's battery fully charged, as the battery operates best in the upper half of its charge. When the battery level drops below 50%, the maximum speed (with a double battery pack) will gradually decrease from 18 km/h to 15 km/h. You can recharge the battery in Scoozy intermittently without causing damage to the battery, unlike an AGM battery.

Read chapters 3.7 and 4 of the manual for additional information on the battery and storing Scoozy.

Why does the range decrease without me driving it?

Just like any mobility scooter, Scoozy has a standby power consumption. In Scoozy, this is visualized in the loss of kilometers without driving. You can charge Scoozy intermittently without damaging the battery but do not leave it continuously connected to the charger. This can eventually draw energy from the battery. Leave the charger connected for a maximum of 12 hours and then remove it.

Can I store Scoozy for an extended period without charging it?

Always charge Scoozy's battery to 100% before storing it for an extended period. Scoozy has special protection against the deep discharge of the battery. If Scoozy's battery is less than 50% charged and Scoozy has not been used for 14 days, the battery goes into a so-called sleep mode. This is displayed as 00 km on the display. You can activate the battery by connecting the charger or by simply driving. After waking up the battery, the correct range with the number of kilometers will be displayed, and you can drive.

We recommend charging the battery every 3 months. When stored for longer periods, the battery may be discharged too deeply, which may require the battery to be replaced.

Why can I only drive 4 km/h when it's under 2°C?

All lithium batteries are damaged when charged at a core temperature below 2 degrees Celsius. Because Scoozy regeneratively brakes, meaning that the battery is charged during braking, Scoozy's speed is limited when the battery is too cold to prevent damage to the battery. A snowflake will appear on the display to alert you to this. Under the menu > info, you can check the battery temperature.

Does Scoozy heat itself up?

Contrary to what is stated in the manual, it turns out in practice that the heating element in Scoozy is not powerful enough to effectively warm up the battery. Unfortunately, there is currently no internal heating available for lithium-ion scooter batteries. Therefore, we advise storing Scoozy in a heated environment, as the battery cools down quickly when stored outside or in an unheated storage area. Watch the video to see how far you can drive with Scoozy at -12 degrees Celsius when starting from a warm place (19 degrees Celsius).

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