

Go Go Carbon ¹	SC25M
Weight Capacity	300 lbs.
Turning Radius	53"
Overall Width ³	17.5"
Overall Length ³	40.3"
Maximum Speed ²	Up to 3.7 mph
Ground Clearance ³	1.5" at motor
Front Wheels	7" solid
Rear Wheels	8" solid
Drivetrain	Rear wheel drive, sealed transaxle, 24-volt brushless
Braking System	Regenerative and electromechanical
Suspension Type	N/A
Seat Sizes	14"W x 11.5"D
Seat-to-Floor Range	21.75"
Per Charge Range ^{3,5}	Up to 10 miles
Standard Electronics	30A Brushless controller
Component Weights ⁶	35 lbs. without battery With battery: 39 lbs.
Battery Requirements ⁴	24V, 12Ah lithium battery
Battery Charger	24V, 2A
Warranty	Frame: lifetime limited Drivetrain: 1-year Electronics: 1-year Battery: 1-year

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 4 AGM or gel cell type required.
- 5 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 6 Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.



PRIDE
MOBILITY

(US) 800-800-8586 • (Canada) 888-570-1113

FOLLOW US      
pridemobility.com

**** Note:** Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to flight to ensure your Go Go Carbon® and its lithium battery are permitted.

GO GO
CARBON

THE FUTURE OF
TRAVEL MOBILITY



FOLLOW US      
pridemobility.com

EXCEPTIONAL PORTABILITY




IT'S ALL IN THE NAME.

Harnessing the power of its full carbon fiber frame, the **Go Go® Carbon** offers exceptional portability in a four-wheel scooter. It weighs just **35 pounds** without its lithium battery, folds down in seconds for convenient transport and storage, and can be towed like luggage through airports and other environments.

- Folds down in just seconds and can be towed like luggage
- Weighs just 35 pounds without the battery
- Full carbon fiber frame, rims, and tiller
- NFC ignition provides extra security
- Vibrant color display on console that displays clear information
- 300-lb. weight capacity
- Top speed of up to 3.7 mph
- Ground clearance of 1.5"
- Range of up to 10 miles per battery charge
- Powered by an airline compliant** lithium battery
- Removable seat back and base for easy cleaning
- Telescoping tiller
- LED lighting for enhanced visibility

COLOR OPTIONS



Black Carbon Fiber

ACCESSORIES

Military Patches



Pride's exclusive lightweight, non-scutting tires

Vibrant full color display

Enhanced LED light